# **Chemistry 20 Solution Stoichiometry - Precipitate Lab**

(3) Complete / neatly done

# (1) **Purpose**:

1 Part B - to produce and collect a precipitate, and then compare theoretical and experimental results

## (5) **Background Information**:

- 1 solution homogeneous mixture of a solute and a solvent
- 1 reactant a chemical used in a reaction
- 1 *product* a chemical produced in a chemical reaction
- 1 precipitate a solid product formed via a chemical reaction
- 1 *filtrate* the solution which goes through the filter paper

# (2) **Experimental Design**:

- 1 manipulated variable amount of sodium carbonate solution
- 1 responding variable amount of precipitate collected

#### (4) **Procedure**:

- 1 Ask the teacher for the required amounts of each solution.
- 1 Rinse pipet before collecting the second solution.
- 1 Weigh filter paper before use.
- 1 Rinse out the reactants beaker into the filter paper.

## (5) **Observations**:

- 1 Mass of filter paper
- 1 Mass of filter paper + precipitate
- 1 Mass of precipitate
- 1 The reaction produced a milky white precipitate.
- Some of the precipitate stuck to the side of the reactants beaker.

#### Conclusion

(5) Part A

$$2 \qquad {\rm K_2CO_{3\,(aq)}} \quad + {\rm CaCl_{2\,(aq)}} \quad \rightarrow \quad 2 \; {\rm KCl_{\,(aq)}} \quad + \qquad {\rm CaCO_{3\,(s)}}$$

3 calculation of theoretical mass

0.50 g of calcium carbonate

# (5) Part B

- 2 Percent error calculation
- 1 Source of error -- some precipitate stuck to the beaker
- 2 Filtrate ions -- chloride and sodium ions

out of 30